

C-Series Rotary Valve



FEATURES

- > 2-year warranty against defects in materials and workmanship
- > Long-life
- > Exceptional Value
- > Available in a variety of inert materials
- > No solenoid valve limitations (leaking, clogging, low CVs)
- > No ceramic valve limitations (seizing)

CONFIGURED OPTIONS

- > Materials
 - Body PEEK
PCTFE
 - Rotary Plug UHMWPE
PTFE
- > Valve Ports
 - 120° 3-port
 - 90° 4-port
 - Distribution 3-way
 - T 3-port
 - T 4-port
 - Loop
 - Distribution 6-way

OEM OPTIONS







- > Custom configurations possible

BASE MODEL

Compatible with C-Series Pumps and Rotary Valve Controllers



SPECIFICATIONS

Description	Materials (Body/Plug)	Catalog Number	Flow Path (Blue)
120° 3-port	PEEK/UHMWPE	10020	
	PEEK/PTFE	10021	
	PCTFE/UHMWPE	10022	
	PCTFE/PTFE	10023	
90° 4-port	PEEK/UHMWPE	10040	
	PEEK/PTFE	10041	
	PCTFE/UHMWPE	10042	
	PCTFE/PTFE	10043	
Distribution 3-way	PCTFE/UHMWPE	10062	
	PCTFE/PTFE	10063	
T 3-port	PEEK/UHMWPE	10115	
T 4-port	PEEK/UHMWPE	10125	
Loop	PCTFE/UHMWPE	10126	
Distribution 6-way	PCTFE/UHMWPE	10106	
	PEEK/UHMWPE	10113	

ROTARY VALVES FOR SYRINGE PUMPS

tricontinent.com



Americas
Tricontinent
Scientific, Inc.

12740 Earhart Ave
Auburn, CA 95602
USA
Tel: +1 800 937 4738 or +1 530 273 8888
liquidhandling.tcs@irco.com

Asia Pacific

Gardner Denver Thomas Pneumatic
Systems (Wuxi), Co., Ltd.
No. 1 New Dong An Road,
Shuofang Town
Wuxi, Xinwu District
Jiangsu 214142 China
Tel: +86 510 6878 2258
tricontinent.cn@irco.com

EMEA

(HEADQUARTERS)
Gardner Denver Thomas GmbH

Livry-Gargan-Str. 10
82256 Fürstenfeldbruck
Germany
Tel: +49 8141 2280 0
thomas.de@irco.com

The information presented in this material is based on technical data and test results of nominal units. It is believed to be accurate and reliable and is offered as an aid to help in the selection of Tricontinent products. It is the responsibility of the user to determine the suitability of the product for the intended use and the user assumes all risk and liability in connection there with. Tricontinent does not warrant, guarantee or assume any obligation or liability in connection with this information.

Photos of products pictured in this catalog do not necessarily represent a specific model number. To obtain further information for custom options, contact your local Tricontinent office.

Printed in USA Form No. MKT90034 E © Tricontinent Scientific, Inc. All rights reserved.

TRICONTINENT
An Ingersoll Rand Business